



Short track to cross

# ACROSS CTO™ **ST**

RX

- High pushability
- Strong 1mm tip
- Enhanced lesion entry
- Low crossing profile
- Distal marker




Tip marker

Short 1 mm tip  
Seamless Tip-to-shoulder

High RBP 17 Atm

Low entry  
profile 0.016"



## Ordering Information / Technical Information

Reference	Balloon Diameter	Balloon Length	Marker	Tip Length	Entry Profile	Max. Guidewire	Usable Length	Min. Guide Catheter / ID
110050483	1,1 mm	5 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"
110100483	1,1 mm	10 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"
110150483	1,1 mm	15 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"
110200483	1,1 mm	20 mm	1	1 mm	0.016"	0.014"	138 cm	5 F / 0.056"

## Control your Size

Balloon Ø	Atm	Atm*	Atm	Atm	Atm	Atm	Atm**	Atm	Atm
1,1 mm	1.04 mm	<b>1.10 mm</b>	1.13 mm	1.16 mm	1.20 mm	1.23 mm	<b>1.28 mm</b>	1.38 mm	1.45 mm

Atm\* Nominal Pressure (NP)      Atm\*\* Rated Burst Pressure (RBP)



The balloon material is a semi-compliant polyamide, providing a 'Work Your Size'.

Short track to cross